



Helm Dental Laboratory

2801 Capital Street
Wylie, TX 75098
972-442-9772

www.HelmDentalLaboratory.com
info@HelmDentalLaboratory.com

ALL information required. Complete all sections. Please keep a copy for your files. We will call you before proceeding if we have any questions.

SCREW-RETAINED DENTURE RX FOR TITANIUM & ACRYLIC PROSTHESIS

Case Information

Arch Selection (Please use separate RX for each arch)

- Maxilla Mandible

Final Screw-retained Denture type

- Titanium & Acrylic

Restorative Doctor's name

Surgical Doctor's Name (for case reference only)

Address

City State Zip

Phone Doctor's e-mail

Patient name Age Male/Female

Doctor's signature

License number / State issued by Today's date

Please send me the following item:

- Rapid Try-In (Appointment #1)
- Wax Try-In (Appointment #2) Shade: _____
- Additional Wax Try-In (Optional Appointment #2A)
- Final Prosthesis (Delivery Appointment)

*See opposite side for summary of appointment steps.
Detailed protocol available upon request.*

Notes:

Please list the implant manufacturer, type & diameter (e.g., Astra Tech EV 4.8, Nobel Active RP 4.3, Straumann BL 4.8, etc.) & tooth number. Additional charges for more than 5 implants.

Implant Mfg.	Type & Diameter	Tooth #

Screw-Retained Prosthesis Package Titanium & Acrylic Ultra includes:


- Rapid Try-In
- Wax Try-In
- Final Screw-Retained Denture
- Model work
- ***Analog not included***
- Bridge screws for final screw-retained denture

Upgrades per arch: (additional charges apply)

- Additional or replacement bridge screws, each
- Restorative chairside assistance, per visit (e.g impressions, try-in, final seating, etc.)
- ThermoguardNightguard
(highly recommended for protection of prosthesis overnight)

Contact Helm Dental Laboratory for RX with current pricing on this product/package.

Please mark missing teeth/ teeth to be extracted with X on tooth chart below.





Summary of Restorative Protocol Screw-Retained Denture Titanium & Acrylic

This protocol outlines the basics of the restorative appointments for completing the delivery of a Titanium & Acrylic Screw-Retained Denture. ***Detailed protocol available upon request.***

RESTORATIVE PHASE

After surgeon releases patient, Restorative doctor sends Screw-Retained Denture RX for Titanium & Acrylic Prosthesis to lab to order a **Rapid Try-In**. The Rapid Try-In is used as a verification index to take a verified master impression.

Allow 10 working days to receive Rapid Try-In.

Appointment 1 - Rapid Try-In (Master Impression and Verification) - Schedule 60 min.

Rapid Try-In appointment is for verification of implant locations.

Items sent by Lab to complete appointment: *Rapid Try-In*

Tools and parts needed to complete appointment: *Tray adhesive, light body impression material, bite registration material, shade tabs, appropriate torque wrench and screw driver.*

Forms: *Screw-Retained Denture RX for Titanium & Acrylic Prosthesis*

- Complete Rapid Try-In process.
- Complete Screw-Retained Denture RX and select Wax Try-In. Indicate shade for final denture.
- Send RX, Rapid Try-In impression, bite, photos and any scans taken to Lab.

Schedule next appointment after 10 working days in lab.

Appointment 2 - Wax Try-In - Schedule 60 min.

Try-in appointment is to verify esthetics and function for final restoration, including final shade.

Items sent by Lab to complete appointment: *Wax Try-In*

Tools & parts needed to complete appointment: *Light body impression material, shade tabs, appropriate torque wrench and screw driver.*

Forms: *Screw-Retained Denture RX, Try-In Checklist, Process & Finish Form.*

- **If there are significant changes** in tooth placement or fit, another Wax Try-In is recommended. This will add an extra appointment to the process. Send RX, Wax Try-In, and Try-In Checklist to Lab.
- **If there are not significant changes**, send RX, Wax Try-In, Try-in Checklist, and signed Process & Finish Form to Lab. Any changes requested on the finished prosthesis after approval of the Try-In will incur additional charges.

Schedule next appointment after 15 working days in lab.

Appointment 3 - Final Prosthesis Delivery - Schedule 30-45 min

Items sent by Lab to complete appointment: *Final Titanium & Acrylic Screw-Retained Prosthesis*

Tools & parts needed to complete appointment: *Light body impression material, Teflon tape, composite, appropriate torque wrench and screw driver.*

Forms: *None*

- Seat **Final Prosthesis**.