



Helm Dental Laboratory

2801 Capital Street

Wylie, TX 75098

972-442-9772

www.HelmDentalLaboratory.com

info@HelmDentalLaboratory.com

ABUTMENT PROVISIONALIZATION PACKAGE ADVANTAGES

- Turns implant dentistry into crown & bridge dentistry
- Provide your referrals with simple & accurate method to take abutment level impression for the final restoration using our exclusive custom pick-up impression coping
- Choice of abutment material and final restoration type

Contact Helm Dental Laboratory for RX with current pricing on this product/package.

NobelProcera® Individualized Abutment, Titanium

- NobelProcera® Titanium abutment, provisional restoration & custom pick-up impression coping

NobelProcera® Individualized Abutment, Gold-shaded

- NobelProcera® Gold-Shaded abutment, provisional restoration & custom pick-up impression coping

NobelProcera® Individualized All Zirconia Abutment

(Zirconia implant connection for anterior teeth only, NP & RP only)

- NobelProcera® Zirconia abutment, provisional restoration & custom pick-up impression coping

OPTIONAL ITEMS

- Custom Seating Guide
- Screwmentable Crown
Provide crown with screw access hole

FINAL RESTORATION TYPE

- Cement-retained
- Screw-retained (check Screwmentable Crown box above if screwmentable provisional is also desired.)

NOTES

- Call me to discuss

All trademarks and company names are the property of their respective owners.

Prescriber and Patient Information - PROVISIONAL RESTORATIONS

Please complete this prescription form, include a copy of the surgical note & send in with your case. Please keep a copy for your files.

Cases will NOT proceed until ALL materials & information are provided.

Doctor's Name

Address

City

State

Zip

Phone

Email

Patient's Name

Age Gender (M/F)

Doctor's Signature

License Number

Today's Date

Final Tooth Shade

IMPLANT INFORMATION

Implant Manufacturer, Type & Diameter	Tooth #	Margin Depths			
		B/F	D	M	L
<i>Lab default margin depths</i>		1.0	0.75	0.75	0.50

Please list the Implant manufacturer, type & diameter (i.e. Nobel Replace™ WP 5.0, NobelActive® RP 4.3, etc.) & tooth numbers of implants. Lab defaults used if no abutment margin depths are indicated.

EMERGENCY PROFILE SELECTION

(Lab default used if no selection made)

- No Tissue Pressure
- Light Tissue Pressure
- Medium Tissue Pressure (Lab default)
- Heavy Tissue Pressure

Bundled Provisionalization package includes:

- Custom designed & milled NobelProcera® Abutment
 - Custom designed & fabricated provisional restoration
 - Custom impression pick-up coping for final impression
 - 1 Abutment screw per implant
 - Nobel Biocare™ Authentication Sticker
 - Model work including soft tissue
- ***Analog not included***

